

Scheduled Patients : Christine Jade Anderson, Able Clinic 3/16/2014 14:47:02

Pos	Status	Patient Name	PID	DOB	Date	Procedure
1	Imaged	Moldenhauer, Joan D	750201	09/07/1942	03/20/2014 15:5	Screening, bilateral, digital, with cad
2	Imaged	Pewe, Pauline C	22563	02/01/1929	03/16/2014 14:4	MAMMOGRAM SCREEN CAD
3	Imaged	Zaghet, Mary N	54304	12/04/1953	03/14/2014 13:5	Screening, bilateral, digital, tomosynt
4	Imaged	Young, Janet E	8624040	12/02/1943	03/14/2014 09:0	Screening, bilateral, digital, with cad
5	Imaged	Waters, Evie D	748127	06/09/1936	03/14/2014 13:4	Screening, bilateral, digital, tomosynt
6	Imaged	Wald, Beverly E / Bev	223235647	09/09/1960	03/14/2014 16:5	Screening, bilateral, digital, with cad
7	Imaged	Wagner, Virginia C / Vik	223249047	07/07/1915	03/14/2014 16:5	Screening, bilateral, digital, tomosynt
8	Imaged	VanHorn, Betty L	223235547	05/05/1960	03/14/2014 13:4	Screening, bilateral, digital, tomosynt
9	Imaged	Ughrin, Kathleen Q / K	51451	09/06/1945	03/14/2014 09:0	Screening, bilateral, digital, with cad
10	Imaged	Tyler, Beverly A / Bev	744617	05/11/1924	03/14/2014 09:0	Screening, bilateral, digital, with cad
11	Imaged	Speulda, Jane I	22409	12/12/1956	03/14/2014 16:5	Screening, bilateral, digital, with cad
12	Imaged	Smith, Theresa P / Ter	223249038	05/15/1919	03/14/2014 13:5	Screening, bilateral, digital, tomosynt
13	Imaged	Sigvaldson, Leslie S	8623466	10/09/1938	03/14/2014 13:4	Screening, bilateral, digital, tomosynt
14	Imaged	Rice-Kramer, Pamela A	3538	02/03/1935	03/14/2014 13:5	Screening, bilateral, digital, tomosynt
15	Imaged	Quinn, Charlotte C / Ch	45396	04/11/1927	03/14/2014 16:5	Screening, bilateral, digital, with cad
16	Imaged	Quigley, Katherine B / K	51312	06/06/1945	03/14/2014 13:5	Screening, bilateral, digital, tomosynt
17	Imaged	Quern, Audrey A	223234950	05/05/1963	03/14/2014 16:5	Screening, bilateral, digital, with cad
18	Imaged	Pewe, Pauline C	22563	02/01/1929	03/14/2014 09:0	Screening, bilateral, digital, with cad
19	Imaged	Peacock, Linda V	752211	10/05/1948	03/14/2014 09:0	Screening, bilateral, digital, with cad
20	Imaged	Nissellus, Beverly M / B	744564	11/11/1924	03/14/2014 09:0	Screening, bilateral, digital, with cad
21	Imaged	Mohr, Beverly A / Bev	44554	07/11/1924	03/14/2014 16:5	Screening, bilateral, digital, with cad

Buttons: Print List, Film Pos, Assign Rad, Show only Imaged, Edit RIS#, Preview Exam Report, Filter By, Sort By, Fetch Images, PreExam, Auto Track, Set Status, Delete, Refresh, Edit Patient, Imaging, Close

PenRad has added a button to allow clients to change an order number, accession number, and/or exam type, for orders canceled through the facility's RIS system after the exam has been imaged or read.

To alter a RIS driven exam, select the exam with a single tap, then tap the Edit RIS button (1) on the reading worklist.

If replacing an order, tap order available on upper list, then tap Assign/Replace on Exam button (2).

Change RIS Procedure Numbers (Accession) for a study : Christine Jade Anderson, 3/16/2014 14:46:39

Pauline C. Pewe DOB: 02/01/1929 AGE: 85 F PID: 22563

Sch.Date	RIS Proc #	Visit #	Procedure	RIS Code	Description	Primary D
03/16/2014 14:...	PT34354		Mammogram	MDB	MAMMOGRAM DIAG BIL	

Manually Add Ris Procedure Order : Christine Jade Anderson, Able Clinic 3/16/2014 14:47:38

Pauline C. Pewe DOB: 02/01/1929 AGE: 85 F PID: 22563

Select Procedure Code

Mammogram	- Mammogram	- MAMMOGRAM SCREEN CAD
MG	- Mammogram	- MAMMO
MDB	- Mammogram	- MAMMOGRAM DIAG BIL
MDR	- Mammogram	- MAMMOGRAM DIAG RIGHT
MDL	- Mammogram	- MAMMOGRAM DIAG LEFT
SCRNB	- Mammogram	- SCREENING MAMMOGRAM
MAMOR	- Mammogram	- DIAGNOSTIC MAMMOGRAM RIGHT
Mammogram	- Mammogram	- MAMMOGRAM SCREEN CAD
MD	- Mammogram	- MAMMOGRAM DIAG
PR0005	- Mammogram	- BREAST TOMOSYNTHESIS STUDY
PR0001	- Mammogram	- DIAGNOSTIC EXAM
PR0002	- Mammogram	- SCREENING EXAM

Pewe, Pauline C: DOB 02/01/1929 : PID 22563 Procedure Filter: Mammogram

RIS Proc # / Accession#: Visit # Date: 03/16/2014 14:47 ? Location Code

Buttons: OK, Cancel, Help

exam/procedure below to assign order to, select RIS order above, then click on "Assign/Replace" button below

Date	Exam Description	RIS Proc Num	RIS Order	Approve Dt	Commit	St
16/2014	Screening, bilateral, digit	PT54234	PT54234			Screening, bilateral, di

Buttons: Show, Show All Locations, Show only available/unused, Print List, Delete, Proc#, Accession:, Search, Refresh, Edit Patient, Prior Exams, Add, Assign/Replace on Exam, Exit, Help

For facilities that create orders within PenRad, if change of order type is required, tap Add button (3) to create new, use drop-down (4) to filter procedure list types to a specified procedure type.

Did you know that PenRad offers a CT lung tracking system also? PenLung allows tracking of these studies, generation of reminder/recall letters, and unresolved exams by radiologist. Offered as an AutoReader system where the results are sent to PenRad via HL7, that displays the narrative report and extracts the key tracking components, and stores the report for reference. Also offered are two budget versions where PenLung is connected to a PACs worklist and stores an exam event when conducted. Then clients are aware that a study requires tracking, eliminating the HL7 result feed, however the report narrative will need to be reviewed and tracked. For low volume sites, we can offer a standalone system. Just enter the patient and track the results.

### Recommendation for product development?

sales@penrad.com | 763.475.3388. Thank you for your business.

© 2017 PenRad Technologies, Inc. 730-23-00261 AB